

## Product Datasheet

# PKA C alpha/beta/gamma (pT198) antibody (orb1187893)

<b>Description</b>	Rabbit polyclonal antibody to PKA C alpha/beta/gamma (pT198)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Bovine, Canine, Human, Mouse, Porcine, Rat, Sheep
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Immunogen</b>	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human PKA C alpha/beta/gamma. The exact sequence is proprietary.
<b>Preservatives</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Concentration</b>	1mg/mL
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P17612</a>
<b>Entrez</b>	<a href="#">5566</a>
<b>Dilution Range</b>	WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500)
<b>Expiration Date</b>	12 months from date of receipt.