

Product Datasheet

Histone H3 (AcK18) antibody (orb1187977)

Description	Rabbit polyclonal antibody to Histone H3 (AcK18)
Species/Host	Rabbit
Reactivity	Bovine, Gallus, Human, Mouse, Porcine, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human Histone H3 (AcK18). The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P68431
Entrez	8350
Dilution Range	WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500)
Expiration Date	12 months from date of receipt.