

## Product Datasheet

### CDK5 (pY15) antibody (orb1187994)

<b>Description</b>	Rabbit polyclonal antibody to CDK5 (pY15)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Bovine, Gallus, Human, Mouse, Porcine, Rat, Sheep
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Immunogen</b>	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human CDK5 (pY15). The exact sequence is proprietary.
<b>Preservatives</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Concentration</b>	1mg/mL
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q00535</a>
<b>Entrez</b>	<a href="#">1020</a>
<b>Dilution Range</b>	WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500)
<b>Expiration Date</b>	12 months from date of receipt.