

Product Datasheet

AMPK gamma 1/2/3 antibody (orb1188257)

Description	Rabbit polyclonal antibody to AMPK gamma 1/2/3
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human AMPK gamma 1/2/3. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P54619
Entrez	5571
Dilution Range	WB (1/500 - 1/1000), IF/IC (1/100 - 1/500)
Expiration Date	12 months from date of receipt.