

Product Datasheet

Glucocorticoid Receptor (pS226) antibody (orb1188263)

Description	Rabbit polyclonal antibody to Glucocorticoid Receptor (pS226)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Glucocorticoid Receptor. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04150
Entrez	2908
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/100)
Expiration Date	12 months from date of receipt.