

Product Datasheet

PRKD1/2/3 (pS738/S742) antibody (orb1188376)

Description	Rabbit polyclonal antibody to PRKD1/2/3 (pS738/S742)
Species/Host	Rabbit
Reactivity	Bovine, Human, Monkey, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human PRKD1/2/3 (pS738/S742). The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q15139
Entrez	5587
Dilution Range	WB (1/500 - 1/1000), IH (1/100 - 1/200)
Expiration Date	12 months from date of receipt.