

Product Datasheet

ATP1A1 (pS16) antibody (orb1188813)

Description	Rabbit polyclonal antibody to ATP1A1 (pS16)
Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human ATP1A1 (pS16). The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P05023
Entrez	476
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.