

## Product Datasheet

### POMGNT1 antibody (orb1189925)

<b>Description</b>	Rabbit polyclonal antibody to POMGNT1
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Recombinant full length protein of human POMGNT1
<b>Preservatives</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Concentration</b>	1mg/mL
<b>Storage</b>	Store at -20°C. Avoid repeated freeze / thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q8WZA1</a>
<b>Entrez</b>	<a href="#">55624</a>
<b>Dilution Range</b>	WB (1/500 - 1/2000), IH (1/50 - 1/200)
<b>Expiration Date</b>	12 months from date of receipt.