

Product Datasheet

REG3 gamma antibody (orb1190125)

Description	Rabbit polyclonal antibody to REG3 gamma
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant full length protein of human REG3 gamma
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q6UW15
Entrez	130120
Dilution Range	WB (1/500 - 1/2000), IF/IC (1/50 - 1/200)
Expiration Date	12 months from date of receipt.