

Product Datasheet

HOXA1/B1/D1 antibody (orb1192893)

Description	Rabbit polyclonal antibody to HOXA1/B1/D1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant full length protein of human HOXA1/B1/D1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P14653
Entrez	3211
Dilution Range	WB (1/500 - 1/2000), IH (1/50 - 1/200)
Expiration Date	12 months from date of receipt.