

Product Datasheet

Mouse IGF1R ELISA Kit (orb1196137)

Description	Mouse Insulin Like Growth Factor 1 Receptor (IGF1R) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q60751
Sensitivity	0.068ng/mL
Entrez	16001
Expiration Date	Please enquire.