

Product Datasheet

Monkey MCP1 ELISA Kit (orb1196223)

Description	Monkey Monocyte Chemotactic Protein 1 (MCP1) ELISA Kit
Reactivity	Monkey
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P61274
Sensitivity	0.054ng/mL
Entrez	102135739
Expiration Date	Please enquire.