

Product Datasheet

Bovine CTSD ELISA Kit (orb1196271)

Description	Bovine Cathepsin D (CTSD) ELISA Kit
Reactivity	Bovine
Range	3.12-200ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	F1MMR6
Sensitivity	0.88ng/mL
Entrez	282883
Expiration Date	Please enquire.