

## Product Datasheet

### General PAF ELISA Kit (orb1196279)

<b>Description</b>	General Platelet Activating Factor (PAF) ELISA Kit
<b>Reactivity</b>	All
<b>Range</b>	0.468-30ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Sensitivity</b>	0.224ng/mL
<b>Expiration Date</b>	Please enquire.