

Product Datasheet

Mouse RSPO3 ELISA Kit (orb1196292)

Description	Mouse R-Spondin 3 (RSPO3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q2TJ95
Sensitivity	0.061ng/mL
Entrez	72780
Expiration Date	Please enquire.