

Product Datasheet

Mouse EGFL7 ELISA Kit (orb1196335)

Description	Mouse EGF Like Domain Protein, Multiple 7 (EGFL7) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9QXT5
Sensitivity	15.5pg/mL
Entrez	353156
Expiration Date	Please enquire.