

## Product Datasheet

### Mouse CD5L ELISA Kit (orb1196341)

<b>Description</b>	Mouse CD5 Antigen Like Protein (CD5L) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.312-20ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9QWK4</a>
<b>Sensitivity</b>	0.136ng/mL
<b>Entrez</b>	<a href="#">11801</a>
<b>Expiration Date</b>	Please enquire.