

Product Datasheet

Human ZG16B ELISA Kit (orb1196379)

Description	Human Zymogen Granule Protein 16 Homolog B (ZG16B) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q96DA0
Sensitivity	0.131ng/mL
Entrez	124220
Expiration Date	Please enquire.