

Product Datasheet

Human WNT5B ELISA Kit (orb1196387)

Description	Human Wingless Type MMTV Integration Site Family, Member 5B (WNT5B) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9H1J7
Sensitivity	0.53ng/mL
Entrez	81029
Expiration Date	Please enquire.