

Product Datasheet

Human WNT11 ELISA Kit (orb1196395)

Description	Human Wingless Type MMTV Integration Site Family, Member 11 (WNT11) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O96014
Sensitivity	0.058ng/mL
Entrez	7481
Expiration Date	Please enquire.