

Product Datasheet

Human VSIG2 ELISA Kit (orb1196400)

Description	Human V-Set And Immunoglobulin Domain Containing Protein 2 (VSIG2) ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q96IQ7
Sensitivity	0.152ng/mL
Entrez	23584
Expiration Date	Please enquire.