

## Product Datasheet

### Human SUZ12 ELISA Kit (orb1196417)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human Suppressor Of Zeste 12 Homolog (SUZ12) ELISA Kit |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 0.312-20ng/mL  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only                                  |
| <b>Application notes</b>   | sandwich   |
| <b>Sample Types</b>        | tissue homogenates or other biological fluids.         |
| <b>Uniprot ID</b>          | <a href="#">Q15022</a>                                 |
| <b>Sensitivity</b>         | 0.113ng/mL   |
| <b>Entrez</b>              | <a href="#">23512</a>                                  |
| <b>Expiration Date</b>     | Please enquire.  |