

Product Datasheet

Human SPLUNC2 ELISA Kit (orb1196424)

Description	Human Short Palate, Lung And Nasal Epithelium Carcinoma Associated Protein 2 (SPLUNC2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q96DR5
Sensitivity	0.059ng/mL
Entrez	140683
Expiration Date	Please enquire.