

Product Datasheet

Human SECTM1 ELISA Kit (orb1196438)

Description	Human Secreted And Transmembrane Protein 1 (SECTM1) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, urine or other biological fluids.
Uniprot ID	Q8WVN6
Sensitivity	0.27ng/mL
Entrez	6398
Expiration Date	Please enquire.