

Product Datasheet

Human PXDN ELISA Kit (orb1196461)

| | |
|----------------------------|--|
| Description | Human Peroxidase Homolog (PXDN) ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Uniprot ID | Q92626 |
| Sensitivity | 0.059ng/mL |
| Entrez | 7837 |
| Expiration Date | Please enquire. |