

Product Datasheet

Human POTE^G ELISA Kit (orb1196467)

Description	Human POTE Ankyrin Domain Family, Member G (POTE ^G) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q6S5H5
Sensitivity	0.117ng/mL
Entrez	404785
Expiration Date	Please enquire.