

Product Datasheet

Human MAP1LC3b ELISA Kit (orb1196518)

Description	Human Microtubule Associated Protein 1 Light Chain 3 Beta (MAP1LC3b) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9GZQ8
Sensitivity	0.068ng/mL
Entrez	81631
Expiration Date	Please enquire.