

Product Datasheet

Human MAP1LC3a ELISA Kit (orb1196519)

Description	Human Microtubule Associated Protein 1 Light Chain 3 Alpha (MAP1LC3a) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9H492
Sensitivity	0.063ng/mL
Entrez	84557
Expiration Date	Please enquire.