

Product Datasheet

Human IFI16 ELISA Kit (orb1196541)

Description	Human Interferon Gamma Inducible Protein 16 (IFI16) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q16666
Sensitivity	0.049ng/mL
Entrez	3428
Expiration Date	Please enquire.