

Product Datasheet

Human CYP2C19 ELISA Kit (orb1196591)

Description	Human Cytochrome P450 2C19 (CYP2C19) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P33261
Sensitivity	0.056ng/mL
Entrez	1557
Expiration Date	Please enquire.