

Product Datasheet

Human CFLAR ELISA Kit (orb1196607)

Description	Human CASP8 And FADD Like Apoptosis Regulator (CFLAR) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O15519
Sensitivity	0.061ng/mL
Entrez	8837
Expiration Date	Please enquire.