

## Product Datasheet

### Human CD27BP ELISA Kit (orb1196616)

<b>Description</b>	Human CD27 Binding Protein (CD27BP) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<b>O15304</b>
<b>Sensitivity</b>	0.068ng/mL
<b>Entrez</b>	<b>10572</b>
<b>Expiration Date</b>	Please enquire.