

Product Datasheet

Rat TRPM6 ELISA Kit (orb1196672)

Description	Rat Transient Receptor Potential Cation Channel Subfamily M, Member 6 (TRPM6) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	F1M7G0
Sensitivity	0.056ng/mL
Entrez	293874
Expiration Date	Please enquire.