

Product Datasheet

Rat PPP2R4 ELISA Kit (orb1196751)

Description	Rat Protein Phosphatase 2A Activator, Regulatory Subunit 4 (PPP2R4) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	B2RYQ2
Sensitivity	0.057ng/mL
Entrez	362102
Expiration Date	Please enquire.