

## Product Datasheet

### Rat PPAR $\gamma$ C1 $\alpha$ ELISA Kit (orb1196754)

<b>Description</b>	Rat Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha (PPAR $\gamma$ C1 $\alpha$ ) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.156-10ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9QYK2</a>
<b>Sensitivity</b>	0.061ng/mL
<b>Entrez</b>	<a href="#">83516</a>
<b>Expiration Date</b>	Please enquire.