

## Product Datasheet

### Rat OXTR ELISA Kit (orb1196773)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rat Oxytocin Receptor (OXTR) ELISA Kit                       |
| <b>Reactivity</b>          | Rat  |
| <b>Range</b>               | 0.156-10ng/mL  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | sandwich   |
| <b>Sample Types</b>        | tissue homogenates, cell lysates or other biological fluids. |
| <b>Uniprot ID</b>          | <a href="#">P70536</a>                                       |
| <b>Sensitivity</b>         | 0.059ng/mL   |
| <b>Entrez</b>              | <a href="#">25342</a>  |
| <b>Expiration Date</b>     | Please enquire.  |