

## Product Datasheet

### Rat MMRN1 ELISA Kit (orb1196815)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rat Multimerin 1 (MMRN1) ELISA Kit             |
| <b>Reactivity</b>          | Rat  |
| <b>Range</b>               | 0.312-20ng/mL                                  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only                          |
| <b>Application notes</b>   | sandwich                                       |
| <b>Sample Types</b>        | tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>          | <a href="#">D4A3E0</a>                         |
| <b>Sensitivity</b>         | 0.113ng/mL                                     |
| <b>Entrez</b>              | <a href="#">500152</a>                         |
| <b>Expiration Date</b>     | Please enquire.                                |