

## Product Datasheet

### Rat GRIA1 ELISA Kit (orb1196872)

<b>Description</b>	Rat Glutamate Receptor, Ionotropic, AMPA 1 (GRIA1) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	15.6-1000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<b>P19490</b>
<b>Sensitivity</b>	6.1pg/mL
<b>Entrez</b>	<b>50592</b>
<b>Expiration Date</b>	Please enquire.