

Product Datasheet

Rat CRHR1 ELISA Kit (orb1196941)

| | |
|----------------------------|--|
| Description | Rat Corticotropin Releasing Hormone Receptor 1 (CRHR1) ELISA Kit |
| Reactivity | Rat |
| Range | 0.312-20ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | tissue homogenates, cell lysates or other biological fluids. |
| Uniprot ID | P35353 |
| Sensitivity | 0.123ng/mL |
| Entrez | 58959 |
| Expiration Date | Please enquire. |