

## Product Datasheet

### Rat AFM ELISA Kit (orb1197015)

<b>Description</b>	Rat Afamin (AFM) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	3.12-200ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<b>P36953</b>
<b>Sensitivity</b>	0.93ng/mL
<b>Entrez</b>	<b>282708</b>
<b>Expiration Date</b>	Please enquire.