

Product Datasheet

Bovine VIP ELISA Kit (orb1197031)

Description	Bovine Vasoactive Intestinal Peptide (VIP) ELISA Kit
Reactivity	Bovine
Range	7.81-500pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P81401
Sensitivity	3.8pg/mL
Entrez	280956
Expiration Date	Please enquire.