

Product Datasheet

Porcine TTR ELISA Kit (orb1197039)

| | |
|----------------------------|--|
| Description | Porcine Transthyretin (TTR) ELISA Kit |
| Reactivity | Porcine |
| Range | 1.56-100ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Uniprot ID | P50390 |
| Sensitivity | 0.68ng/mL |
| Entrez | 397419 |
| Expiration Date | Please enquire. |