

Product Datasheet

Porcine TNNI3 ELISA Kit (orb1197046)

Description	Porcine Troponin I Type 3, Cardiac (TNNI3) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	A5X497
Sensitivity	12.6pg/mL
Entrez	100049696
Expiration Date	Please enquire.