

Product Datasheet

Porcine SELE ELISA Kit (orb1197068)

Description	Porcine Selectin, Endothelium (SELE) ELISA Kit
Reactivity	Porcine
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P98110
Sensitivity	0.237ng/mL
Entrez	397508
Expiration Date	Please enquire.