

## Product Datasheet

### Porcine RANTES ELISA Kit (orb1197079)

<b>Description</b>	Porcine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	3.12-200ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q2EN89</a>
<b>Sensitivity</b>	1.15ng/mL
<b>Entrez</b>	<a href="#">396613</a>
<b>Expiration Date</b>	Please enquire.