

## Product Datasheet

### Bovine RANTES ELISA Kit (orb1197080)

<b>Description</b>	Bovine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	78.1-5000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">O97919</a>
<b>Sensitivity</b>	18.8pg/mL
<b>Entrez</b>	<a href="#">327712</a>
<b>Expiration Date</b>	Please enquire.