

## Product Datasheet

### Porcine MPO ELISA Kit (orb1197107)

<b>Description</b>	Porcine Myeloperoxidase (MPO) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	0.781-50ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">K7GRV6</a>
<b>Sensitivity</b>	0.35ng/mL
<b>Entrez</b>	<a href="#">100517120</a>
<b>Expiration Date</b>	Please enquire.