

Product Datasheet

Porcine MMP9 ELISA Kit (orb1197109)

Description	Porcine Matrix Metalloproteinase 9 (MMP9) ELISA Kit
Reactivity	Porcine
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q2VI03
Sensitivity	0.56ng/mL
Entrez	654325
Expiration Date	Please enquire.