

Product Datasheet

Bovine MIP1b ELISA Kit (orb1197123)

Description	Bovine Macrophage Inflammatory Protein 1 Beta (MIP1b) ELISA Kit
Reactivity	Bovine
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q17QA1
Sensitivity	9.8pg/mL
Entrez	414347
Expiration Date	Please enquire.