

Product Datasheet

Bovine MIP1a ELISA Kit (orb1197124)

Description	Bovine Macrophage Inflammatory Protein 1 Alpha (MIP1a) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8SQA6
Sensitivity	7.4pg/mL
Entrez	282170
Expiration Date	Please enquire.